

# CERTIFICATE OF ANALYSIS

For

**Lithium Buffer 3**

Part Number: 80-2099-83

Batch (Control) Number: W27915

Quantity: 2 Litres

Expiry Date: 13 March 2028

The above product gave the following results on analysis:

pH = 3.148 @ 20.0°C

pH Specification = 3.13 - 3.17 @ 20.0°C

Conductivity = 34.3 mScm<sup>-1</sup> @ 20.0°C

Conductivity Specification = 33.4 - 39.7 mScm<sup>-1</sup> @ 20.0°C

Chemical QC Signature:

Date of Issue: 19 March 2025

Neil Booth

|                                                                                                                                       |                    |                    |                        |
|---------------------------------------------------------------------------------------------------------------------------------------|--------------------|--------------------|------------------------|
| Postal Address                                                                                                                        | Telephone          | Telefax            | Email:                 |
| Biochrom Ltd<br>Unit 7, Enterprise Zone<br>3970 Cambridge Research Park<br>Beach Drive, Waterbeach<br>Cambridge<br>CB25 9PE, England. | +44 (0)1223 423723 | +44 (0)1223 420164 | support@biochrom.co.uk |
| Web Site: <a href="http://www.biochrom.co.uk">http://www.biochrom.co.uk</a>                                                           |                    |                    |                        |

Registered in England No: 3526954 Registered Office: Unit 7, Enterprise Zone, 3970 Cambridge Research Park, Beach Drive, Waterbeach, Cambridge, CB25 9PE England.